The Nuclear Fuel Cycle, Lab

In the first part of the laboratory activities, the assistant will demonstrate different detectors for measuring ionizing radiation. In the second part the students are supposed to work more or less without supervision. Prepare by considering how the task shall best be performed.

The following text describes how the laboratory experiments are to be performed. Read in the textbook about how the different detectors work in general.

- (8.2.3 Pulse shape and dead time 202-204
 - 8.3 Gas counters 204-206
- 8.3.3 GM-counters 210-211
 - **8.4** Semi-conductor detector 212-213, 215-217
 - **8.5** Scintillation detectors 218-222
- 8.7.5 Multi-channel analyzers 226
- **8.7.6** *γ*-spectroscopy 226-228

Ion chamber / Geiger-Müller counter

Three different Geiger-Müller counters will be demonstrated

- Determine the β_{max} of a sample by first determining the half-thickness of an absorber (see appendix)
- > Determine the measuring efficiency of the detector for β- and γ-radiation. Assume that the detector is working in the ion chamber region (even though it's not)
- > Remember to take the background into account.

Scintillation detector

Determine the activity [Bq] of a sample. A solution with known specific activity will be provided for the isotope.

Multi-channel analyzer

- ▶ Record the spectrum of Eu-152
- Determine the measuring efficiency of the detector at three different energies, using data for Eu-152 provided in the appendix.

Mats Jansson e-mail: matsja@kth.se phone: 790 7193



Fig. 6 Halveringstjockleken för betapartiklar som funktion av energin

Commis	sariat à l'Energie Atomique 1920 - 1935 - 293 - 293 - 293 - 293 - 293 - 293 - 293 - 293 - 293 - 293 - 293 - 293	Bureau National de Métrologie
	Référence Multigamma ¹⁵² Eu	
ource étalon type EGMA 2 ate de référence : le 05.09.	№ 5179 à 12 h 00(1)	
ENERGIE ⁽²⁾ KeV	Taux d'émission photonique (3) ۲/همو/4π	Incertitude (4) %
121,78 ± 0,02	1,168.10 ⁴	2,4
244,69 ± 0,02 244,7 ⁽⁵⁾ 251,76 ± 0,08	3,072.10 ³	2,6
295,94 ± 0,03	1,751.10 ²	3,7
344,28 ± 0,02	1,093.10 ⁴	2,1
367,76 ± 0,03	3,479.10 ²	4,2
$\begin{array}{r} 411,12 \pm 0.02 \\ 411,3 {}^{(5)} \\ 416,04 \pm 0.07 \end{array}$	9,524.10 ²	3,1
$\begin{array}{r} 443,98 \ \pm \ 0.02 \\ 444,0 \ \pm \ 0.10 \end{array} \qquad $	1,273.10 ³	2,2
$\begin{array}{r} 564.02 \pm 0.04 \\ 566.36 \pm 0.08 \end{array} $	2,538.10 ²	6;1.
686,4 — 688,69 ± 0,08	3,516.10 ²	3,5
778,90 ± 0,03	5,383.10 ³	2,3
867,38 ± 0,03	1,723.10 ³	2,0
$963,43 \pm 0,08 964,05 \pm 0,03 964,0 (5)$	5,996.10 ³	2,2
$1085,83 \pm 0,03 \\ 1086,4 ^{(5)} \\ 1089,72 \pm 0,04$	4,903.10 ³	2,7
1109,2 — 1112,0 (5)	5,611.10 ³	2,1

6,048.10²

7,218.102

8,571.103

1298,7 (5)

I.P. Nº 2 - 91190 - GIF-SUR-YVETTE - FRANCE

1112,08 \pm 0,04

1212,94 ± 0,04

1292,86 ± 0,07

1299,13 ± 0,04

1408,03 ± 0,03

Tél. 941.80.00 - Télex 69641 F Energat Saciay

2,8

3,6

2,0



Appendix